

Biomark HD

- **Time savings.** Perform 9,216 reactions in as little as 30 minutes.
- **Cost savings.** Do more experiments at 10x less cost per datapoint.
- **Flexibility.** Use commercially available assays and DNA binding dyes.

The Biomark™ HD system sets a new standard for high-throughput real-time PCR, end-point PCR and digital PCR, with benefits that are impossible to reproduce using many other conventional PCR systems. Integrated fluidic circuit (IFC) technology both prepares and performs thousands of reactions in nanoliter volumes, saving time and money and reducing pipetting steps by 95%. The Fluidigm® Biomark HD, together with IFCs and the IFC Controller for loading samples and assays, streamlines workflows for applications demanding sensitivity and dynamic range at an extremely high throughput, such as single-cell analysis.



HIGH-THROUGHPUT DETECTION



The system integrates thermal cycling and detection of PCR assays for all Dynamic Array™ IFCs and Digital Array™ IFCs. It acquires data for each reaction chamber on the IFC simultaneously and can operate in either real-time or end-point detection mode for gene expression or for genotyping experiments, respectively.

ASSAYS

The Biomark HD offers an open and flexible platform accommodating reagents and chemistries of choice. Also, the entire system, from the footprint of the IFCs to the architecture of analysis and database software, adheres to industry standards and ensures integration with established workflows.

ANALYSIS SOFTWARE

The Biomark HD is bundled with data collection and data analysis software. Real-Time PCR Analysis software displays the analyzed data in multiple formats, including color-coded maps of every reaction chamber on the IFC, amplification curves, and numeric tables. Results may be easily managed, annotated and archived.



Real-Time PCR Analysis software generates the heat map and amplification plot.

GENE EXPRESSION ASSAYS


Delta Gene™ assays for mid-plex gene expression studies provide robust, high-quality real-time PCR gene expression assays. The assays enable users to take full advantage of the Biomark HD with minimal experimental setup time using validated protocols that provide quality results.

- Amplicons designed to cross an intron whenever possible to avoid genomic DNA amplification
- Designed to any RefSeq including human, mouse, rat, microbe, and plant (minimum of 48 assays per order)
- Single-cell gene expression protocols available
- Turnaround time of three weeks for bioinformatically tested assays and six weeks for wet-lab tested assays
- Custom panels/pathways designed upon request

GENOTYPING ASSAYS


SNP Type™ assays provide a high-throughput, low-cost single-nucleotide polymorphism (SNP) genotyping solution which enables rapid assay design and polymorphism screening. The assays are based on allele-specific PCR SNP detection chemistry and combine the advantages of minimum experimental setup time and flexible assay choice with the reliability of Dynamic Array IFCs.

- Designed to target species with available sequence information
- Three to four week design and turnaround time with customer-provided sequences (minimum of 24 assays per order)
- Access to loci-specific primer sequences assures reproducibility
- Compatible with specific target amplification (STA) protocol for improving results from samples of low quality and/or concentration, or from species with large genome sizes (>human); necessary STA primers provided




Call map view for 48 cattle samples and 48 SNP Type™ assays (left); cluster plot for a typical SNP Type assay (right)

C_q difference between 1, 10, and 100 cells easily distinguished



Single T_m Peak



1, 10, and 100 cells, custom EvaGreen® Assay linearity data (triplicates shown)

The Biomark HD is compatible with multiple system components to meet a variety of application and sample throughput needs.

IFC CONTROLLER COMPATIBILITY

	IFC Controller MX	IFC Controller HX	IFC Controller RX
Gene expression	48.48 Dynamic Array IFC	96.96 Dynamic Array IFC	—
SNP genotyping	48.48 Dynamic Array IFC	96.96 Dynamic Array IFC	192.24 Dynamic Array IFC
Digital PCR	12.765 Digital Array IFC 48.770 Digital Array IFC qdPCR 37K IFC	—	—
Experiment tracking	Barcode		
Gas pressure	Internal compressor		
Interface	USB and ethernet		
IFC Controller MX, HX, RX or WX software	Touchscreen interface for operating and tracking		
Dimensions (approx.)	19 x 9.5 x 13 inches; 48.5 x 24 x 33 cm		

IFCs FOR GENOTYPING

	48.48 Dynamic Array IFC	96.96 Dynamic Array IFC	192.24 Dynamic Array IFC	FLEXsix IFC
Assay transfer rate	98.00%			
Call rate	99.00% or greater			
Call accuracy	99.75% or greater			
Dimensions	SBS compatible (128 mm x 85 mm x 14 mm)			
Inlet spacing on input frame	4.5 mm pitch			
Liquid transfer steps	96	192	216	Variable
Assay inlets	48	96	24	6x12
Sample inlets	48	96	192	6x12
Reaction chambers	2,304	9,216	4,608	864
Instrument compatibility	EP1™ reader, IFC Controller MX, FC1™ reader, Biomark HD	EP1, IFC Controller HX, FC1, Biomark HD	EP1, IFC Controller RX, FC1, Biomark HD	EP1, IFC Controller HX, FC1, Biomark HD

IFCs FOR GENE EXPRESSION

	48.48 Dynamic Array IFC	96.96 Dynamic Array IFC
Quantitative resolution	Two-fold difference in starting copy with 99.7% confidence and 6-log of dynamic range	
Inlet spacing on input frame	4.5 mm pitch	
Dimensions	SBS compatible (128 mm x 85 mm x 14 mm)	
Liquid transfer steps	96	192
Assay inlets	48	96
Sample inlets	48	96
Reaction chambers	2,304	9,216
Reaction volume	9 nL	6.7 nL
Instrument compatibility	Biomark HD, IFC Controller MX	Biomark HD, IFC Controller HX

IFCs FOR DIGITAL PCR

	12.765 Digital Array IFC	48.770 Digital Array IFC qdPCR 37k IFC
Detection sensitivity	Single copy (if copy is present in the reaction chamber)	
Dimensions	SBS compatible (128 mm x 85 mm x 14 mm)	
Inlet spacing on input frame	4.5 mm pitch	
Minimum input volume/sample	8 µL (12 samples per array)	4 µL (48 samples per array)
Liquid transfer steps	12	48
Sample inlets	12	48
Reactions per sample	765	770
Total reaction chambers	9,180	36,960
Individual reaction volume	6 nL	0.85 nL
Total reaction volume/sample	4.6 µL (per sample)	0.65 µL (per sample)
Instrument compatibility	Biomark HD, EP1 reader, IFC Controller MX	

SOFTWARE SPECIFICATIONS

The Biomark software suite was designed to offer a simple and intuitive user interface while continuing to offer all key data analysis features required by today's scientists. To simplify and expedite data analysis, the software suite includes key features:

- **Heat Maps**
IFC layout maps have color-coded reactions to identify C_i or delta C_i trends or SNP genotyping calls.
- **Multiple IFC Scoring**
Scoring combines multiple experiments to decrease analysis time and improve accuracy for rare genotypes.
- **Sample and Assay Mapping**
Sample and assay information can be pasted directly from Microsoft® Excel® and imported from a database of a previously saved template.
- **Quality Scores**
Quality scoring allows quick and efficient sorting through PCR curves and cluster calls to identify and exclude those that do not meet criteria.

SYSTEM COMPONENTS	
Excitation filters (center-width, in nm)	485-20, 530-20, 580-25 (two empty positions)
Emission filters (center-width, in nm)	525-25, 570-30, 645-75 (two empty positions)
Thermal Control	4–99 °C range Heating (65–90 °C) >2 °C/sec Cooling (90–65 °C) >1 °C/sec
Software	Fluidigm Real-Time PCR Analysis Fluidigm Genotyping Analysis Fluidigm Digital PCR Analysis Fluidigm Data Collection

Corporate headquarters

7000 Shoreline Court, Suite 100
South San Francisco, CA 94080 USA
Toll-free: 866 359 4354 | Fax: 650 871 7152
fluidigm.com

Sales

North America | +1 650 266 6170 | info-us@fluidigm.com
Europe/EMEA | +33 (0)1 60 92 42 40 | info-europe@fluidigm.com
Japan | +81 (0)3 3662 2150 | info-japan@fluidigm.com
China (excluding Hong Kong) | +86 (0)21 3255 8368 | info-china@fluidigm.com
Asia | +1 650 266 6000 | info-asia@fluidigm.com
Latin America | +1 650 266 6000 | info-latinamerica@fluidigm.com

©2014 Fluidigm Corporation. All rights reserved. Fluidigm, the Fluidigm logo, Biomark, Dynamic Array, Digital Array, Delta Gene, SNP Type, EP1, and FC1 are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

For Research Use Only. Not for use in diagnostic procedures.

PN 100-4158 B1